YDENUPET^{SP}

DEUTSCH

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YDENUPET^{SP}

ENGLISH

FOR THE DENUDATION AND PIPETTING OF OOCYTES AND EMBRYOS

INDICATIONS FOR USE

YDENUPET SP is used for the removal of cumulus oocyte complexes (COC) from an oocyte prior to Intracytoplasmic Sperm Injection (ICSI) and *In Vitro* Fertilization (IVF) as well as for the handling of oocytes and embryos during assisted reproductive techniques (ART). YDENUPET SP is not intended for biopsy of cells from oocytes or embryos.

PRODUCT DESCRIPTION

The tip of YDENUPET SP is available in a range of internal diameter sizes. Sizes 125 - 175 µm are suitable for oocyte denudation, sizes 200 - 600 µm are for oocyte and embryo handling procedures.

PRODUCT SPECIFICATIONS

- Material: Polycarbonate, non-pyrogenic
- Primary packaging: Polypropylene
- Secondary packaging: Cellulose
- Total volumetric capacity (depending on the size): 20 – 25 µL
- Proximal dimension: OD 900 µm
- Endotoxins: < 0.05 EU/device
- Sterility: sterile (SAL 10⁻⁶)
- MEA 1-Cell to expanded blastocyst in 96h: ≥80%
- A certificate of analysis and MSDS are available upon request.

STERILE R Sterilized by gamma irradiation

CONTRAINDICATION

No known contraindication.

WARNINGS

- Only for intended use.
- Do not use if package is damaged.
- Not sterile after opening.
- Do not re-sterilize.
- Do not reuse.
- Multiple use may lead to in cross-contamination.
- Each laboratory should decide for itself which YDENUPET SP diameter to use for each specific procedure! Internal diameter of the tip must be larger than the diameter of the sample. If the appropriate tip is not selected, sample damage is possible during ART procedures.
- YDENUPET SP should be used only with pipettors compatible with the proximal dimensions of YDENUPET SP (0.9 mm).

PRECAUTIONS

- VITROMED recommends handling under aseptic conditions.
- The products may only be disposed of in accordance with the applicable regulations for waste products.

Rx Only Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

YDENUPET SP should only be used for its intended purpose by staff trained with the procedure.

General handling:

- ① Attach a YDENUPET SP carefully on the pipetting device.
- ② Rinse the YDENUPET SP by the pipetting device.
- ③ Before working with the oocytes, embryos or blastocysts, aspirate a little medium into the YDENUPET SP. Only then start with the removal of the COC or the handling of oocytes and embryos.
- ④ To avoid any air bubbles, do not completely empty fluid from the YDENUPET SP when using it.
- ⑤ Dispose of the YDENUPET SP after use by appropriate measures.

Denudation:

- ① Transfer the oocytes into media with an enzyme constituent, hyaluronidase, which promotes the denuding process.
- ② Aspirate/expel repeatedly through pipette tips of incrementally smaller internal diameters to remove the outer layers of cumulus and corona cells.
- ③ Always keep the media level visible in the YDENUPET SP.

STORAGE INSTRUCTIONS AND STABILITY

- Do not use after expiry date.
- Store at 15 – 25 °C.
- Keep away from sunlight and heat.
- Keep dry.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com. Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established or treated.

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VITROMED GmbH