# HYADASE

## DEUTSCH

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## ENGLISH

**HYADASE** is a ready-to-use media for the denudation of oocytes prior to the intracytoplasmic sperm injection (ICSI). It contains 80 IU/mL of Hyaluronidase in HEPES buffered Human Tubal Fluid (HTF) Medium with Human Serum Albumin (HSA). Hyaluronidase catalyses the enzymatic cleavage of hexosamine bonds between hyaluronic acid and glucosaminoglycans, which anchor cumulus cells to the oocyte, and thereby facilitates removal of cumulus cells during the denudation process.

### COMPOSITION

- HEPES HTF
- Gentamicin
- Hyaluronidase 80 IU/mL
- Human Serum Albumin

### PRODUCT SPECIFICATIONS

- Endotoxins: < 1.00 EU/mL
- Sterility: sterile (SAL 10<sup>-3</sup>)
- pH: 7.20 - 7.40
- Osmolality: 257 - 273 mOsm/kg
- Mouse Embryo Assay (Blastocysts after 96 h): ≥ 80%
- A certificate of analysis and MSDS are available upon request.

**STERILE A** Sterilized using aseptic processing techniques.

### CONTRAINDICATION

No known contraindication.

### PRE-USE CHECK

- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discoloured, turbid or a precipitate is observed.
- Do not use if seal of the container is damaged upon delivery.

### INSTRUCTIONS FOR USE

**HYADASE** should be used only for its intended purpose by qualified and trained persons.

- 1) Prepare a petri dish containing one drop **HYADASE** and several drops of **HEPES** (80 µL to 100 µL each) and cover with **VOIL**. Equilibrate the dish to 37 °C.
- 2) Transfer 3 to 6 cumulus oocyte complexes (COC) into the **HYADASE** droplet and gently pipette the COCs up and down using a denuding pipette (**DENUPET**) to remove the oocytes from the cumulus cells.
- 3) After a maximum of 20 to 30 seconds transfer the partially denuded oocytes with a fine denuding pipette (**DENUPET**) into a **HEPES** droplets and remove the remaining granulosa cells by further denudation.
- 4) Wash the oocytes thoroughly in the remaining **HEPES** droplets.
- 5) Oocytes are now ready for the ICSI procedure.

### STORAGE INSTRUCTIONS AND STABILITY

- Store **HYADASE** between 2 °C to 8 °C.
- Do not freeze before use.
- Leave product in its original packaging.
- Do not use after expiry date.
- Keep away from sunlight.

### WARNING AND PRECAUTIONS

**VITROMED** recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5). Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed colour. Discard minimal amounts of excess medium remaining in the bottle. Handle all samples as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling samples. This product contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus

(HIV) by approved testing methods. Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, the product carries an extremely remote risk for transmission of viral diseases. The theoretical risk for the transmission of CJD is also considered extremely remote. No cases of viral disease or CJD have ever been identified for albumin. **HYADASE** contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction). Only for intended use. Ensure correct disposal.

### BOVINE SOURCED HYALURONIDASE

The pharmaceutical grade hyaluronidase used in this product is derived from bovine testis and is certified with a TSE risk evaluation Certificate of Suitability (CEP). The animals from which the hyaluronidase is derived, are determined "fit for human consumption" and originate from countries with "negligible BSE risk", as determined Resolution "Recognition of the Bovine Spongiform Encephalopathy Risk Status of Member Countries", adopted by the OIE. According to the WHO guidelines on Tissue Infectivity Distribution in Transmissible Spongiform Encephalopathies (2010), testes from bovine source are classified as Category IC material (i. e. "Tissues with no detected infectivity").

### Rx Only

Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician.

## DE / EN

# HYADASE

Stück pro VE / Piece per unit	REF
5 x 1 mL Röhrcen / vial	V-HYA

Technischer Support  
Technical support

**VITROMED GmbH**  
Raiffeisenstr. 15a  
40764 Langenfeld  
Germany

+ 49 2173 . 200 41 30  
support@vitromed.com  
www.vitromed.com

Stand der Informationen 01-Aug-2019  
Date of Information