БЪЛГАРСКИ
ČEŠTINA
CPICKИ
DEUTSCH
EΛΛHNIKA
ENGLISH
ESPAÑOL
FRANÇAIS
HRVATSKI
ITALIANO
MAGYAR
NEDERLANDS
POLSKI
PORTUGUÊS
ROMÂNĂ
SHQIP
SLOVAK
SLOVENŠČINA

VITROMED GmbH
Hans-Knöll-Str. 6
07745 Jena
Germany
+49 36 41 . 5 39 19 76
support@vitromed.com
www.vitromed.com

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
ΓΛΩΣΣΑΡΙΟ ΣΥΜΒΟΛΩΝ (ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements)

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FOR THE EXAMINATION OF SPECIMEN

INDICATIONS FOR USE

IVF Labware is intended for the examination of specimen, e.g. sperm, oocytes and tissue, derived from the human body.

- **VCENTERWELL:** For preparing, examining the cells and observing the growth of the cells under the same conditions.
- **VCULTUREDISH:** For the cultivation of cells, searching for Cumulus-oocyte-complexes in follicular fluid and observation of cell growth under the microscope.
- **VICSIDISH:** For examining the quality of a sample under the microscope, e.g. sperm, enables a better image under the microscope at a higher magnification.
- **V4WELLDISH:** For the handling of gametes and observation of the cell growth.
- **VCRYOVIAL:** For cryopreservation of cells and tissues for *in vitro* diagnostic purposes.
- **VROUNDTUBE:** For andrological application, examination, and diagnosis of infertility.
- **VCONICALTUBE:** For centrifuging of samples and examining them after centrifugation.
- **VSAMPLE CONTAINER:** For storing sperm samples for further *in vitro* diagnostic procedures.
- **VSEROLOGICAL PIPET:** For aspirating and dispensing liquids with manual and electronic pipette controls.
- **VTRANSFER PIPET:** For manual aspiration and dispensing of liquids.
- **V TIP:** For the transmission of liquids.

PRODUCT DESCRIPTION

IVF Labware are ready-to-use plastic disposables.

PRODUCT SPECIFICATIONS

- **Material:**
 - **VCENTERWELL:** Polystyrene.
 - **VCULTUREDISH:** Polystyrene.
 - **VCRYOVIAL:** Polypropylene.
 - **VICSIDISH:** Polystyrene.
 - **V4WELLDISH:** Polystyrene.
 - **VROUNDTUBE:** Polystyrene tube with LDPE cap.
 - **VCONICALTUBE:** Polystyrene tube with LDPE cap.
 - **VSAMPLE CONTAINER:** Polypropylene.
 - **VSEROLOGICAL PIPET:** Polystyrene.
 - **VTRANSFER PIPET:** LDPE.
 - **V TIP:** Polypropylene.
- **Sterility:** sterile (SAL 10⁻⁴)
- **Endotoxins:** < 20 EU/device
- **MEA 1-Cell to expanded blastocyst in 96h:** ≥80%
- **Human Sperm Survival Assay (HSSA):** ≥ 70 %
- A certificate of analysis and MSDS are available upon request.

 Sterilized using irradiation.

CONTRAINDICATION

No known contraindication.

WARNINGS

- Only for intended use.
- Do not use after expiry date.
- Do not use if package is damaged.
- Not sterile after opening.
- Do not re-sterilize.
- Do not reuse.
- Multiple use may lead to in cross-contamination.

PRECAUTIONS

- VITROMED recommends handling under aseptic conditions.
- The products may only be disposed of in accordance with the applicable regulations for waste products.

INSTRUCTIONS FOR USE

IVF Labware should only be used for its intended purpose by staff trained with the procedure.

VCENTERWELL

1. Fill the center well and the outer chamber with equilibrated media according to your protocol. For the center well compartment, the recommended volume is 1 mL when used with oil, otherwise 2 mL. For the outer well compartment, the recommended volume is 2 mL when used with oil, otherwise 4 mL.
2. Dispose after use by appropriate measures.

VCULTUREDISH:

1. For small volume cell culturing use VCULTUREDISH-35, for large volume cell culturing VCULTUREDISH-90, in case of droplet culture place small droplets (20-50 µL) on the bottom of the dish.
2. Overlay the drops with oil.
3. VITROMED recommends preincubating the medium and Oil when working with CO₂ dependent media.
4. Place the cells into the culture medium.
5. Incubate the cells according to your protocol
6. Dispose after use by appropriate measures.

VICSIDISH:

1. Place small droplets on the bottom of the dish. Recommended volume per drop is 5 – 50 µL depending on the intended use.
2. Overlay the drops with oil. The recommended volume is 5 mL.
3. Place the cells to examine into the droplet and examine them e.g. with a microscope.
4. Dispose after use by appropriate measures.

V4WELLDISH:

1. To culture the cells, fill the wells with medium and place the dish in the incubator or heating unit before use.
2. The recommended volume of medium per well is 0,5 mL in combination with oil, otherwise 1.2 mL. VITROMED recommends the use of oil.
3. Dispose after use by appropriate measures.

VCRYOVIAL:

1. Handle the semen according to your Semen Freeze protocol.
2. The max. volume of the freezing suspension that can be placed into the tube is 2 mL.
3. Dispose after use by appropriate measures.

VROUNDTUBE:

1. Fill the sample into the tube.
2. Loosely put the cap over the tube for gas changing.
3. Dispose after use by appropriate measures.

VCONICALTUBE:

1. Fill the prepared sample into the tube.
2. Close the cap tightly before centrifugation. The maximum centrifugation speed is 4800 g (15 minutes).
3. Dispose after use by appropriate measures.

VSAMPLE CONTAINER:

1. Fill the sample in the container.
2. Screw the cap properly to avoid leaking.
3. Dispose after use by appropriate measures.

VSEROLOGICAL PIPET:

1. Attach a VSEROLOGICAL PIPET on the pipetting device.
2. Place the tip of the pipette into the preferred liquid and aspirate it using the pipetting device.
3. Follow the scale on the pipette.
4. The pipette can be used for mixing chemical solutions, mixing cell suspension, or transferring liquids. When required volume is drawn up, gently tap the side of the pipette to remove excess liquid.
5. Hold pipette at 10 – 20° angle when dispensing the liquid. Do not blow through the pipette to remove the liquid.
6. Dispose after use by appropriate measures.

VTRANSFER PIPET:

1. Gently squeeze the bulb of the pipette, then place the pipette upright in the liquid.
2. Gently release the bulb to draw the liquid into the pipette. When the required volume has been drawn up, gently tap the side of the pipette to remove excess liquid.
3. Hold the pipette at a 10 – 20° angle and squeeze the bulb to dispense the liquid.
4. Dispose after use by appropriate measures.

V TIP:

1. Attach a V TIP on the pipetting device.
2. Following the instructions of your pipetting device.
3. Hold the pipette vertically and dip the tip only 2-3 mm into the liquid.
4. Dispose after use by appropriate measures.


STORAGE INSTRUCTIONS AND STABILITY

- Store at room temperature.
- Keep away from sunlight and heat.
- Keep dry.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com. Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established.

SYMBOL GLOSSARY (ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements)

Symbol (reference number)	Title of symbol
 (5.1.1)	Manufacturer
 (5.1.4)	Use-by date
 (5.1.5)	Batch code
 (5.1.6)	Catalog number
 (5.2.4)	Sterilized using irradiation
 (5.2.6)	Do not re-sterilize
 (5.2.8)	Do not use if package is damaged
 (5.2.11)	Single sterile barrier system
 (5.2.14)	Single sterile barrier system with protective packaging outside
 (5.3.1)	Fragile, handle with care
 (5.3.2)	Keep away from sunlight
 (5.4.2)	Do not reuse
 (5.4.3)	Consult instructions for use
 (5.5.1)	<i>In Vitro</i> Diagnostic Medical Device

Directive 98/79/EC on *in vitro* diagnostic medical devices

IVF Labware	
Piece per unit	REF
V-CENTERWELL, 300 x 1	V-CWD-60
V-CULTUREDISH-60, 300 x 1	V-CTD-60
V-CULTUREDISH-90, 300 x 1	V-CTD-90
V-CULTUREDISH-35, 300 x 1	V-CTD-35
V-ICSIDISH-50, 300 x 1	V-ICD-50
V-4WELLDISH, 200 x 1	V-4WD
V-CRYOVIAL-2, 20 x 50	V-CRV-2
V-ROUNDTUBE-14, 500 x 1	V-RTU-14
V-ROUNDTUBE-5, 500 x 1	V-RTU-5
V-CONICALTUBE-15, 500 x 1	V-CTU-15
V-SAMPLE CONTAINER-100, 300 x 1	V-SAC-100
V-SEROLOGICAL PIPET-1, 300 x 1	V-SRO-1
V-SEROLOGICAL PIPET-5, 200 x 1	V-SRO-5
V-SEROLOGICAL PIPET-10, 200 x 1	V-SRO-10
V-TRANSFER PIPET-3, 500 x 1	V-TRP-3
V-TIP 1000, 500 x 1	V-TIP-1000
V-TIP 200, 500 x 1	V-TIP-200

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