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ENGLISH

FOR SEPARATION OF MOTILE SEMEN FROM SPERM SAMPLES

PRODUCT DESCRIPTION AND INTENDED USE

°GRAD is a sterile colloidal suspension of silica particles stabilized with covalently bound hydrophilic silanes. Available as 100% stock solution and diluted in HEPES buffered Human Tubal Fluid (HTF) medium as 40% or 80% solution. °GRAD100 is a stock solution for preparing a density gradient system for semen preparation. °GRAD40 and °GRAD80 are gradient systems for semen preparation.

COMPOSITION

• Silane-coated Silica • HEPES HTF
Gentamicin and/or Phenol Red can be added upon request.

PRODUCT SPECIFICATIONS

- Endotoxins: < 0.5 EU/mL • Sterility: sterile (SAL 10⁻³)
- pH: 7.20 - 7.90 (Release criteria: 7.20 - 7.60)
- Density: 1.097 - 1.107 g/mL (°GRAD 80 %)
1.115 - 1.125 g/mL (°GRAD 100 %)
< 1.65 cP (°GRAD 80 %)
< 1.75 cP (°GRAD 100 %)
- Viscosity: 310 - 340 mOsm/kg (°GRAD 40 %)
310 - 340 mOsm/kg (°GRAD 80 %)
300 - 330 mOsm/kg (°GRAD 100 %)
- Human Sperm Survival Assay tested ≥ 80% survival after 4 h exposure of spermatozoa to the test medium.
- A certificate of analysis, MSDS and SSCP are available upon request.

STERILE A Sterilized using aseptic processing techniques.

CONTRAINDICATION

No known contraindication.

PRE-USE CHECK

- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discoloured, turbid or a precipitate is observed.
- Do not use if seal of the container is damaged upon delivery.
- Do not use if the expiry date has been exceeded.

INSTRUCTIONS FOR USE

°GRAD should be used only for its intended purpose by qualified and trained persons.

PREPARATION OF GRADIENTS

VITROMED advises to prepare a dual gradient system (45% - 90% or 40% - 80%) starting from °GRAD 100%. If preferred a multi-layer gradient can also be used (e. g. 45% - 70% - 90%).

- ① Mix °GRAD 100% thoroughly before use.
- ② Prepare a 90% gradient medium by mixing 1 part HEPES buffered HTF medium (e. g. °HEPES PLUS) to 9 parts °GRAD 100%.
- ③ Prepare a 45% gradient medium by mixing 5.5 parts HEPES buffered HTF medium (e. g. °HEPES PLUS) to 4.5 parts °GRAD 100%.

Note: Alternatively, any HEPES-buffered medium that is recommended for sperm preparation can be used for the preparation of the gradient media.

- ④ Mix well to get homogeneous solutions.
- ⑤ Follow the instructions „DENSITY GRADIENTS FOR SPERM SELECTION“ in order to use the prepared density gradient media for sperm selection.

Note: Gradient media should be prepared under sterile conditions (e. g. LAF bench ISO Class 5). For optimal results, prepare the gradient media a maximum of 24 hours prior to use, store at 2 - 8 °C and warm gradients to room temperature or 37 °C one hour before use.

DENSITY GRADIENTS FOR SPERM SELECTION

- ① Mix °GRAD bottles by 5 bottle inversions before use.
- ② Warm up all components and samples to room temperature or 37 °C.
- ③ Transfer 2.5 mL of °GRAD 80% into a sterile disposable centrifuge tube.
- ④ Carefully layer 2.5 mL of °GRAD 40% on top of the °GRAD 80% layer, without disturbing the phase interface. Pay attention that the two layers stay distinctly separated.
- ⑤ Gently place up to 2.5 mL of liquefied semen on top of the °GRAD 40% layer.
- ⑥ Centrifuge for 15 - 18 minutes at between 350 and 400 x g. Warning: Should a pellet not be visible after the initial centrifugation step, perform a second centrifugation step for 3 to 5 minutes.
- ⑦ Remove the supernatant leaving approximately 0.3 mL of fluid above the pellet.
- ⑧ Add 2 - 3 mL of sperm washing medium (e. g. °SPERM WASH), resuspend the pellet and transfer to a new centrifuge tube.
- ⑨ Centrifuge for 8 - 10 minutes at 300 x g.
- ⑩ Remove the supernatant and repeat steps 8 and 9.
- ⑪ Remove the supernatant and resuspend the pellet in the desired volume of sperm washing medium (e. g. °SPERM WASH).
- ⑫ The sample is now ready for use.

Note: For highly viscous samples centrifuge at up to 500 x g. Take care not to exceed this centrifugal force.

STORAGE INSTRUCTIONS AND STABILITY

- Store V^{GRAD} between 2 °C to 8 °C.
- Do not freeze before use.
- Leave product in its original packaging.
- V^{GRAD} should be used for a period of up to 7 days after opening, provided that it is handled under aseptic conditions and the products are stored at 2-8°C.
- The shelf life of the product is maintained after transport (max. 5 days) at higher temperatures (≤ 37 °C).
- Keep away from sunlight.

WARNING AND PRECAUTIONS

Not sterile after opening. Do not re-sterilize. Do not reuse. Multiple use may lead to cross-contamination. VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5). Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed colour. Discard minimal amounts of excess medium remaining in the bottle. Handle all samples as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling samples. V^{GRAD} contains gentamicin (depending on the product variant), appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction). Only for intended use. Disposal may only be carried out in accordance with the regulations for waste products.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com. Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the competent authority of the Member State in which the user and/or patient is established.

SYMBOL GLOSSARY

(ISO 15223-1:2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements)

SYMBOL (reference number)	TITLE OF SYMBOL
 (5.1.1)	Hersteller Manufacturer
 (5.1.4)	Verfallsdatum Use-by date
 (5.1.5)	Chargennummer Batch code
 (5.1.4)	Katalognummer Catalog number
 (5.2.2)	Sterilisiert mittels aseptischer Verfahrenstechniken Sterilized using aseptic processing techniques
 (5.2.8)	Nicht verwenden, wenn Packung beschädigt Do not use if package is damaged
 (5.3.2)	Vor Sonnenstrahlen schützen Keep away from sunlight
 (5.3.7)	Bei 2 – 8 °C lagern Store at 2 – 8 °C temperature limits
 (5.4.2)	Nicht zur Wiederverwendung Do not reuse
 (5.4.3)	Gebrauchsanweisung beachten Consult instructions for use
	Medizinprodukt Medical Device
	Eindeutige Produktkennung Unique Device Identifier

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Technischer Support **VITROMED GmbH**