V DENU PET HANDLE

DEUTSCH

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ENGLISH

V DENU PET HANDLE is used in combination with denuding pipettes during assisted reproductive treatment for removal of cumulus cells from cumulus-oocyte complexes (COC), handling of oocytes, pronuclear stages, embryos and blastocysts.

MATERIALS INCLUDED

- V DENU PET HANDLE
- Grip (6): Optionally, a grip can be attached for better handling. To do this, slide the enclosed grip (6) over the blue shaft (1) of the V DENU PET HANDLE.
- Spare parts (O-Rings): To replace the O-rings, unscrew and remove the clamp sleeve (2). Remove the two black O-rings and the white spacer. Place the new black O-ring with the small hole against the barrel, then the white spacer, and finally the O-ring with the larger hole. Replace the clamp sleeve (2) and tighten it.
- Allen key

RECOMMENDATION

Before using the V DENU PET HANDLE for the first time, it is recommended to practice on discarded tissue (e. g. cumulus masses).

INSTRUCTIONS FOR USE

V DENU PET HANDLE should be used only for its intended purpose by qualified and trained persons.

Attaching a denuding pipette

- 1 Loosen the clamp sleeve (2) and depress the plunger (4) until stainless-steel plunger (3) protrudes 0.5 to 1.0 cm above the clamp sleeve (2).
- 2 Gently attach the denuding pipette to the stainless-steel plunger (3) until it stops against the "O" rings inside the barrel.
- 3 Tighten the clamp sleeve (2).
- 4 To rinse the pipette, press the plunger (4) until it reaches the stop, immerse the pipette into a drop of media and slowly release the plunger (4). Expel the media by depressing the plunger (4) as before.

Setting the pipetting volume

- 1 The piston skirt is preset to 3 µL. To reduce the delivery volume, turn the piston shaft clockwise until the desired sample volume is opposite the reference mark on the spring housing.
- 2 A full turn reduces the capacity by 0.2 µL.

Denuding/handling

- 1 To draw-up cells, depress the plunger (4), immerse the denuding pipettes in the culture medium and slowly release the plunger (4).
- 2 To release the cells, depress the piston until it reaches the stop.

CLEANING

- Clean the surface of the V DENU PET HANDLE by wiping with isopropyl alcohol or another approved cleaning agent.
- V DENU PET HANDLE cannot be autoclaved.

STORAGE

Store the V DENU PET HANDLE in a clean, dry place at room temperature.


WARNING AND PRECAUTIONS

- Only for intended use.
- V DENU PET HANDLE is compatible with denudation pipettes with a proximal dimension of 900 µm outer diameter and 600 µm inner diameter.

DE / EN

V DENU PET HANDLE

Stück pro VE / Piece per unit	REF
1	V-DH02

VITROMED
Germany

Technischer Support
Technical support
Assistance technique
Assistenza tecnica
Asistencia técnica

VITROMED GmbH
Hans-Knöll-Str. 6
07745 Jena
Germany

+49 36 41 . 5 39 19 76

support@vitromed.com
www.vitromed.com

Stand der Informationen
Date of information

14.01.2022