DEUTSCH

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ENGLISH

FOR THE COLLECTION, HANDLING, MICROMANIPULATION AND TRANSFER OF OOCYTES AND EMBRYOS

INDICATIONS FOR USE

*HEPES PLUS is intended for oocyte collection, intracytoplasmic sperm injection (ICSI) fertilization, and oocyte and embryo washing and handling procedures. *HEPES PLUS is also intended for transfer of embryos to the uterine cavity.

CONTRAINDICATION

No known contraindications.

PRODUCT DESCRIPTION

*HEPES PLUS is a ready-to-use HEPES buffered Human Tubal Fluid (HTF) Medium.

COMPOSITION

- HEPES HTF
- Gentamicin
- Bicarbonate
- Phenol Red
- Lactate
- Ultra pure water
- Amino acids
- Glucose
- Human Serum Albumin
- Physiological Salts

PRODUCT SPECIFICATIONS

- Endotoxin (LAL) (USP<85>): < 0,25 EU/mL
- Sterility (USP <71> / EP 2.6.1.): sterile (no growth)
- pH: 7,00 - 7,40
- MEA 1-cell to expanded blastocyst in 96 h after 120-minute exposure: ≥80%
- A certificate of analysis, MSDS and SSCP are available upon request.

STERILE A Sterilized using aseptic processing techniques.

WARNING

- Do not use after expiry date.
- Do not re-sterilize after opening.
- VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5).
- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discolored, turbid or a precipitate is observed.
- Do not use if seal of the bottle, or the bottle (sterile barrier) is damaged upon delivery.
- Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed color.
- *HEPES PLUS contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).
- *HEPES PLUS is not intended for follicle flushing procedures.

PRECAUTIONS

- Prescription only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Only for intended use.
- Do not freeze before use.
- Leave product in its original packaging.
- Discard minimal amounts of excess medium remaining in the bottle.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.
- *HEPES PLUS contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods. Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, the product carries an extremely remote risk for transmission of viral diseases. The theoretical risk for the transmission of CJD is also considered extremely remote. No cases of viral disease or CJD have ever been identified for albumin.
- The product may only be disposed of in accordance with the applicable regulations for waste products.

INSTRUCTIONS FOR USE

*HEPES PLUS should be used only for its intended purpose by qualified and trained persons. VITROMED recommends that each laboratory uses its own optimized and validated laboratory procedures.

Oocyte Collection

- ① Warm ³HEPES PLUS to 37 °C prior to use. Do not incubate in an incubator with increased CO₂ concentration.
- ② Remove the desired volume of product, sufficient to flush the equipment and for the procedure, using aseptic techniques. VITROMED recommends daily fresh aliquotation in sterile tubes with well closing lids.
- ③ Place ³HEPES PLUS in sterile collection tubes to receive retrieved oocytes.
- ④ Use your own optimized and validated laboratory procedures for oocyte collection and follicle aspiration.
- ⑤ After collection, the oocytes can be transferred to a suitable culture medium for further IVF procedures or used for micromanipulation.

Microinjection (ICSI)

- ① Warm ³HEPES PLUS to 37 °C in culture dishes or flasks prior to use. Do not incubate in an incubator with increased CO₂ concentration.
- ② Before microinjection place collected oocytes to remove cumulus cells into the wells or droplets on a dish loaded with hyaluronidase. After Denudation incubate the oocytes in culture medium until microinjection.
- ③ For the microinjection, transfer the denuded oocytes into an oil-coated dish with microdroplets from ³HEPES PLUS.
- ④ Fertilize the oocyte by ICSI with prepared sperm and transfer the oocyte to an appropriate culture medium (e. g. ³ONESTEP).

Embryo Transfer

³HEPES PLUS can be used for embryo transfer. VITROMED recommends not to change the medium, and to use culture medium containing the embryo for transfer.

STORAGE INSTRUCTIONS AND STABILITY

- Do not use after expiry date.
- Store between 2 °C to 8 °C.
- The product can be used for a period of up to 7 days after opening, provided that it is handled under aseptic conditions.

FEEDBACK AND OBLIGATION TO INFORM

Thank you for purchasing our product. If you have any suggestions for us to improve our products, please send them to qc@vitromed.com. Report any serious incidents that have occurred with the product immediately by email to qc@vitromed.com and the responsible competent authority.

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Technischer Support
Technical support

VITROMED GmbH
Raiffeisenstr. 15a
40764 Langenfeld
Germany

+ 49 2173 . 200 41 30
support@vitromed.com
www.vitromed.com

Stand der Informationen / Date of information 12-Sep-2022