V-PENTOXIFYLLINE

DEUTSCH

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ENGLISH

V-PENTOXIFYLLINE is a HEPES buffered Human Tubal Fluid (HTF) Medium with Pentoxifyllin. It is intended to distinguish dead from non-moving sperm by inducing a short-term motility, in testicular tissue (TESE) or necrozoospermic ejaculates. For use in research procedures only.

COMPOSITION

- HEPES HTF
- Sodium Hyaluronate
- Human Serum Albumin
- Gentamicin
- Physiological Salts
- Taurine
- Pentoxifylline

PRODUCT SPECIFICATIONS

- Product specifications, a certificate of analysis and MSDS are available upon request.

STERILE Sterilized using aseptic processing techniques

CONTRAINDICATION

No known contraindication.

PRE-USE CHECK

- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discoloured, turbid or a precipitate is observed.
- Do not use if seal of the container is damaged upon delivery.

INSTRUCTIONS FOR USE

V-PENTOXIFYLLINE should be used only for its intended purpose by qualified and trained persons.

1. Warm V-PENTOXIFYLLINE to 37 °C, without CO₂ equilibration.
2. Examine the sperm sample in a mixture consisting of one part sample and one part V-PENTOXIFYLLINE (For example: Add 20µl V-PENTOXIFYLLINE to 20µl sample).
3. Incubate for 10 min. The effect of V-PENTOXIFYLLINE last for 60 min.
4. Sperm evaluation should occur on a heated microscope stage, at 37°C.

STORAGE INSTRUCTIONS AND STABILITY

- Store V-PENTOXIFYLLINE between 2 °C to 8 °C.
- Do not freeze before use.
- Leave product in its original packaging.
- Do not use after expiry date.
- Keep away from sunlight.

WARNING AND PRECAUTIONS

VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5). Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed colour. Discard minimal amounts of excess medium remaining in the bottle. Handle all samples as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling samples. This product contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods.

Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, the product carries an extremely remote risk for transmission of viral diseases. The theoretical risk for the transmission of CJD is also considered extremely remote. No cases of viral disease or CJD have ever been identified for albumin.

V-PENTOXIFYLLINE contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction). Only for intended use. Ensure correct disposal.

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Technischer Support
Technical support

VITROMED GmbH
Raiffeisenstr. 15a
40764 Langenfeld
Germany

+49-2173-2004130
support@vitromed.com
www.vitromed.com

Stand der Informationen
Date of information

01-Aug-2020

IFU-V-PENTOXIFYLLINE_V01