## DEUTSCH

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## ENGLISH

<sup>Y</sup>PVP is a ready-to-use 7% or 10% pharma grade polyvinylpyrrolidone (PVP) solution containing HEPES buffered Human Tubal Fluid (HTF) medium and Human Serum Albumin (HSA). <sup>Y</sup>PVP is used for reducing the motility and to simplify the immobilization and ease of handling sperm during intracytoplasmic sperm injection (ICSI). <sup>Y</sup>PVP also prevents sperm sticking to the injection pipette (<sup>Y</sup>ICSIPIPET) and gives better control of media and sperm in the ICSI pipette.

### COMPOSITION

- Polyvinylpyrrolidone
- Pyruvate
- HEPES HTF
- Human Serum Albumin
- Lactate
- Glucose
- Physiologicall Salts
- Gentamicin

### PRODUCT SPECIFICATIONS

- Endotoxins: < 0.5 EU/ml
- Sterility: sterile (SAL 10<sup>-3</sup>)
- pH: 7.25 - 7.45
- Mouse Embryo Assay (Blastocysts after 96 h): ≥ 80%
- A certificate of analysis and MSDS are available upon request.

**STERILE A** Sterilized using aseptic processing techniques.

### CONTRAINDICATION

No known contraindication.

### PRE-USE CHECK

- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if discoloured, turbid or a precipitate is observed.
- Do not use if seal of the container is damaged upon delivery.

### INSTRUCTIONS FOR USE

<sup>Y</sup>PVP should only be used for its intended purpose by qualified and trained persons.

- 1 Place several droplets of HEPES buffered media (5 µL - 10 µL) and a droplet of <sup>Y</sup>PVP (5 µL - 10 µL) in a petri dish, cover with <sup>Y</sup>OIL and warm up to 37 °C.
- 2 Add 1 µL - 2 µL of washed sperm to the centre of the <sup>Y</sup>PVP droplet.
- 3 Immobilize and aspirate the sperm for oocyte injection during ICSI.

### STORAGE INSTRUCTIONS AND STABILITY

- Store <sup>Y</sup>PVP between 2 °C to 8 °C.
- Do not freeze before use.
- Leave product in its original packaging.
- Do not use after expiry date.
- Keep away from sunlight.

### WARNING AND PRECAUTIONS

VITROMED recommends only handling of media with aseptic techniques (e. g. LAF-bench ISO Class 5). Do not use medium that shows any evidence of particulate matter, is cloudy, or has changed colour. Discard minimal amounts of excess medium remaining in the bottle. Handle all samples as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling samples. This product contains albumin, a derivative of human blood. All donors used in its manufacture were individually tested and found to be nonreactive for hepatitis B surface antigen (HBsAg) and antibodies to hepatitis C virus (HCV) and human immunodeficiency virus (HIV) by approved testing methods. Donors of the source material have been screened for Creutzfeldt-Jakob disease (CJD). Based on effective donor screening and product manufacturing processes, the product carries an extremely remote risk for transmission of viral diseases. The theoretical risk for the transmission of CJD is also considered extremely remote. No cases of viral disease or CJD have ever been identified for albumin. <sup>Y</sup>PVP contains gentamicin, appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction). Only for intended use. Ensure correct disposal.

**Rx Only** Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician.

## DE / EN

<b><sup>Y</sup>PVP</b>	
Stück pro VE / Piece per unit	REF
5 x 0,5 mL Röhrrchen / 5 x 0,5 mL vial	V-PVP-7
5 x 0,5 mL Röhrrchen / 5 x 0,5 mL vial	V-PVP-10

Technischer Support  
Technical support

**VITROMED GmbH**  
Raiffeisenstr. 15a  
40764 Langenfeld  
Germany

+ 49 2173 . 200 41 30  
support@vitromed.com  
www.vitromed.com

Stand der Informationen 18-Oct-2019  
Date of information