

The logo consists of three blue circles arranged horizontally, each containing a white letter: 'I', 'V', and 'F'.

IVF

LABWARE

A VITROMED BRAND

A large, clear petri dish is shown in the upper left corner, partially cut off by the edge of the frame.


Labware Catalogue

Our IVF Labware brand offers
secure solutions during your
laboratory procedures

www.ivflabware.com


A circular logo containing the CE mark, indicating compliance with European standards.

CE





IVF LABWARE of VITROMED is a specially designed portfolio of high quality medical plastic ware and particularly **IVF** tested for IVD use.



IVF
LABWARE
A VITROMED BRAND

Focus

on Quality Control

All our IVF Labware products are:

- **Cleanroom** manufactured
- **Individually soft blister packed**
- **Particle & scratch free**
- **Non toxic**
- **Endotoxin** tested
- 96h 1-cell Mouse Embryo Assay (**MEA**) > 80% tested
- **HSSA** > 70% tested
- Polymers with **USP Class 6** standard
- **Bisphenol A** free
- **Sterile by radiation**

V-CENTERWELL

Center-well-dish

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- Single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Center well Petridish	V-CWD-60	60 x 15	300



individual soft blister packed

V-CULTUREDISH-35, 60, 90

Petridishes with lid in different sizes

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Petridish 35 mm	V-CTD-35	35 x 10	300
Petridish 60 mm	V-CTD-60	60 x 15	300
Petridish 90 mm	V-CTD-90	90 x 15	300



individual soft blister packed

V-ICSIDISH-50

50 mm petridish

Details

- Clear virgin polystyrene with excellent optical properties
- 300 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
ICSI Petridish	V-ICD-50	50 x 9	300



individual soft blister packed

V-4 WELLDISH


4-well-dish

Details

- 4 well dish with flat bottom, volume 1.8 mL
- Wells marked with numbers
- Clear virgin polystyrene with excellent optical properties
- 200 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
4-Welldish	V-4WD	66 x 66	200



individual soft blister packed

V-ROUNDTUBE-5

Round bottom tube 5 mL

Details

- Printed graduations
- Dual position snap-cap
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box,
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Round bottom tube 5 ml	V-RTU-5	75 x 12	500



individual soft blister packed

V-ROUNDTUBE-14

Round bottom tube 14 mL

Details

- Printed graduations
- Dual position snap-cap
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Round bottom tube 14 ml	V-RTU-14	100 x 17	500



individual soft blister packed

V-CONICALTUBE-15

15 mL tube with conical shape

Details

- Printed graduations with white label area
- The practical screw cap prevents the liquid from leaking
- Clear virgin polystyrene with excellent optical properties
- 500 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Centrifuge tube	V-CTU-15	120 x 17	500



individual soft blister packed

V-CRYOVIAL-2

2 mL vials for cryopreservation with screw cap

Details

- Round bottom, screw cap with white silicon ring
- Bottom specially designed for use in cryogenic racks and trays, free standing
- Polypropylene tubes to withstand temperatures to -196°C
- Printed graduations with white label area
- 50 pcs/bag, 1000 pcs/box

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/Unit
Cryovial 1.8 ml	V-CRV-2	47 x 12	1000



V-SAMPLE CONTAINER

100 mL lockable sample container

Details

- Blue screw cap
- Graduations up to 90 mL, total volume 100 mL
- Clear virgin polypropylene with excellent optical properties
- 300 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Sperm samplecup	V-SAC-100	65 x 58	300



individual soft blister packed

V-SEROLOGICAL PIPET-1, 5, 10

Serological pipettes of different sizes

Details

- Printed +/- graduations
- color code for size identification
- Clear virgin polystyrene with excellent optical properties
- 300 & 200 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Serological pipet 1 ml	V-SRO-1	280 x 5	300
Serological pipet 5 ml	V-SRO-5	340 x 8	200
Serological pipet 10 ml	V-SRO-10	340 x 10	200



individual soft blister packed

V-TRANSFER PIPET-3

3 mL pasteur (transfer) pipette

Details

- Small tip ensures consistent drop size
- total volume of 3 mL
- Low density polyethylene
- 500 pcs/box
- single soft blister packed

Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
Pasteur pipete 3 ml	V-TRP-3	150 length	500



individual soft blister packed

V-TIP-200, 1000

High quality pipette tips in different sizes

Details

- color code for size identification: 200 µl - yellow, 1000 µl - blue
- 500 pcs/box
- single soft blister packed


Ordering Information

Type	Ref No.	Dimension (mm)	Qty/unit
200 µl pipet tip	V-TIP-200	n.a.	500
1000 µl pipet tip	V-TIP-1000	n.a.	500



individual soft blister packed





IVF LABWARE

A VITROMED BRAND



Our IVF LABWARE manufacturing facility



JENA

Hans-Knöll-Str. 6
07745 Jena
Germany
Phone: +49 3641 539 19 76



other VITROMED
manufacturing facility:
LANGENFELD, Germany

other VITROMED offices:
US, Greece, Spain



VITROMED
Germany

www.ivflabware.com
order@vitromed.com

Connect with **VITROMEDonline**



For additional product or technical information please contact us at support@vitromed.com