



# VITROMED

Germany

## Product Catalogue Europe

We stand for German precision,  
high quality products and  
compassionate customer service

[www.vitromed.com](http://www.vitromed.com)







## WHO WE ARE

VITROMED was founded in 2011 and is now a leading global company in manufacturing high-quality medical devices used during all stages in fertility treatment in IVF clinics. Our advanced medical device manufacturing facilities are located in Germany, Jena and Langenfeld, but our products are available worldwide. In 2020 and 2021 VITROMED established representative offices in Spain and USA for marketing and other service operations.

### Our Commitment to Quality & Service

VITROMED is an ISO 13485 certified company that ensures the manufacturing of top-notch products through a quality management system and with the usage of advanced technologies and the latest research. Most VITROMED products are MEA and LAL tested and undergo stringent, lot-wise quality control before product release.

Most of VITROMED products are CE marked and FDA registered. More information can be found on our website.

### Ask our Scientific Advisors

We have well-trained and highly professional technical and service teams, who are always ready to resolve any problem our customers face and provide information regarding our products.

We strive to supersede expectations and take extreme care in meeting customer requirements. Over the years, we have been providing support with authentic products and services to numerous clinics around the globe. We ensure our customer satisfaction with high-quality customized products augmented with beautiful design and integrity.

See what our clients say about us at [www.vitromed.com](http://www.vitromed.com).

Additionally, we have a team of experts, headed by our Scientific Advisors, the late Dr. Patrick Quinn and Dr. Ralf Böhm, who are leading several research and development projects.



VITROMEDIA


MICROTOOLS

OTHER PRODUCTS





# VITROMEDia





## V-OIL

V-OIL is a ready to use light paraffin oil for covering tissue culture in *in vitro* fertilization, embryo culture and micro-manipulative procedures.

### Details

- Endotoxin tested ( $< 0.1$  EU/mL)
- 120h 1-cell time lapse morphokinetics mouse embryo assay (MEA)  $\geq 80\%$
- Blastocyst cell count on day 6 of mouse embryo development
- Peroxide value (POV) determination ( $< 0.1$  mEq/KG)
- Prewashed twice with high-purity water
- Available in 100 mL and 500 mL with PETG and Glass bottles

### Ordering Information

Type	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
Paraffin Oil	V-OIL-P100	100	PETG	1
Paraffin Oil	V-OIL-P500	500	PETG	1
Paraffin Oil	V-OIL-G100	100	GLASS	1
Paraffin Oil	V-OIL-G500	500	GLASS	1



# MICROTOOLS



CE



### V-ICSIPET

V-ICSIPET HOLDING and INJECTION pipettes are manufactured from special cleaned borosilicate glass under cleanroom conditions.

#### Details

- Endotoxin tested (< 0.05 EU/device)
- 96h 1-cell Mouse Embryo Assay (MEA) ≥ 80%
- Gamma sterilized
- Customized single packaging for easy and secure removal
- Color-coded according to purpose
- 1 x 10 pipettes individually, sterile packed

#### Ordering Information \*

Type	Ref Nr.	Dimension (µm)	Angle (°)	Qty/Unit
INJECTION	V-INJ-S3-30	7 (OD)	30	10
INJECTION	V-INJ-S3-35	7 (OD)	35	10
HOLDING	V-HOLD-30	120 (OD)	30	10
HOLDING	V-HOLD-35	120 (OD)	35	10

\* For custom dimensions please contact our support team



### V-DENUPET

V-DENUPET pipettes are manufactured from high quality medical polymer under cleanroom conditions.

#### Details

- Endotoxin tested (< 0.05 EU/device)
- 96h 1-cell Mouse Embryo Assay (MEA) ≥ 80%
- Gamma sterilized
- Easy and secure removal
- Unique color code for each size
- V-DENUPET - 5 x 10 pipettes

#### Ordering Information

Type	Ref Nr.	Dimension (µm)	Qty/Unit
DENUPET	V-DEN-	125, 135, 140, 150, 175 200, 275, 300, 600	50



### V-ICSIPET POLAR BODY BIOPSY

POLAR BODY BIOPSY pipettes are manufactured from special cleaned borosilicate glass under cleanroom conditions.

#### Details

- Endotoxin tested (< 0.05 EU/device)
- 96h 1-cell Mouse Embryo Assay (MEA) ≥ 80%
- Gamma sterilized
- Customized single packaging for easy and secure removal
- Color-coded according to purpose
- 1 x 10 pipettes individually, sterile packed

#### Ordering Information \*

Type	Ref Nr.	Dimension (µm)	Angle (°)	Qty/Unit
Polar Body	V-PGDP-V1-30	23 (OD)	30	10
Polar Body	V-PGDP-V1-35	23 (OD)	35	10

\* For custom dimensions please contact our support team

#### Info

➤ Please check the product availability in your region while ordering.





# OTHER PRODUCTS





**V-DENUPET**  
HANDLE


V-DENUPET HANDLE is used in combination with denuding pipettes during assisted reproductive treatment for removal of cumulus cells from cumulus-oocyte complexes (COC), handling of oocytes, pronuclear stages, embryos and blastocysts.

**Details**

- V-DENUPET HANDLE allows secure volume adjustments till 3 µL
- Precise control for manipulation of oocytes and embryos
- Smooth and responsive plunger making it ideal for daily lab routine
- Easy and quick fixation of denudation pipettes with different sizes

**Ordering Information**

Type	Ref Nr.	Material	Qty/Unit
Handset	V-DH02	Metal	1



**V-SPERM**  
MORPHO SLIDE

V-SPERM MORPHO SLIDE provide a pure and clear exposure of sperm for qualitative assessment of morphology for *in vitro* diagnostic purposes in IVF/ART.

**Details**

- Stains: cresyl violet acetate and new methyl blue
- Time and cost efficient
- Ready to use
- 1 x 10 pre-stained slides and 10 cover slides

**Ordering Information**

Type	Ref Nr.	Stains	Qty/Unit
Pre-stained Slides	V-SMS-01	Cresyl violet acetate & new methylene blue	10 (Slides and Cover slides)



**V-GLASSPET**

V-GLASSPET pipettes are designed for easy liquid handling.

**Details**

- Endotoxin tested (< 0.05 EU/device)
- 96h 1-cell Mouse Embryo Assay (MEA) ≥ 80%
- Gamma sterilized
- Plugged with sterile cotton or unplugged, sizes: 150 mm and 230 mm
- 2 mL volume
- 10 x 25 pipettes, sterile blister packed

**Ordering Information**

Type	Ref Nr.	Dimension (mm)	Classification	Qty/Unit
Pasteur Pipette	V-GPT-150	150 (Length)	Plugged	250
Pasteur Pipette	V-GPT-150NP	150 (Length)	Non-Plugged	250
Pasteur Pipette	V-GPT-230	230 (Length)	Plugged	250
Pasteur Pipette	V-GPT-230NP	230 (Length)	Non-Plugged	250





**V-MARKER**  
PEN

V-MARKER PEN with its ultra fine tip are ideal for the marking of dishes, cryovials and straws.

**Details**

- Endotoxin tested (< 0.05 EU/device)
- 96h 1-cell Mouse Embryo Assay (MEA) ≥ 80%
- Non-toxic and xylene free
- Available in black, red and blue
- 1 x 4 pens in your choice of colors

**Ordering Information**

Type	Ref Nr.	Classification	Qty/Unit
Marker Pen	V-PEN-DC	Red, Blue & Black	4
Marker Pen	V-PEN-BK	Black	4
Marker Pen	V-PEN-BL	Blue	4
Marker Pen	V-PEN-R	Red	4



**V-CARE**  
STERI CONTROL

V-CARE STERI CONTROL enables an easy control of the aseptic working environment and medium handling. The product enables an unmatched level of control over the sterile liquid handling processes used in ART. V-CARE SteriControl ensures the traceability of sterile working conditions.

**Details**

- Tryptic Soy Broth (TSB) medium
- Growth promotion tested
- Shelf life 12 months from date of production

**Ordering Information**

Type	Ref Nr.	Volume (mL)	Packaging	Qty/Unit
STERI CONTROL	V-CARE-SC20P	20	PETG	1
STERI CONTROL	V-CARE-SC20G	20	GLASS	1

CE

IVF  
LABWARE  
A VITROMED BRAND



Visit [www.ivflabware.com](http://www.ivflabware.com)  
and our **IVF LABWARE** product catalogue for more details





# OUR DISTRIBUTORS



- our loyal contracted partners of trust -



Scan me

## Our corporate offices



**JENA**

Hans-Knöll-Str. 6  
07745 Jena  
Germany  
Phone: +49 3641 539 19 76



**LANGENFELD**

Raiffeisenstr. 15a  
40764 Langenfeld  
Germany  
Phone: +49 2173 200 41 30

Connect with **VITROMEDonline**



For distributor listing, please visit us at [www.vitromed.com/distributors](http://www.vitromed.com/distributors)  
For additional product or technical information please contact us at [support@vitromed.com](mailto:support@vitromed.com)