

CERTIFICATE OF ANALYSIS

LOT NUMBER: 103VM3602
PRODUCT NAME: V-HYADASE
REFERENCE: V-HYA
EXPIRY DATE: 15-Mar-2022

R_x Only Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician.

TEST	SPECIFICATION	RESULT
pH	7.20-7.40	7,35
OSMOLALITY	257 - 273 mOsm/kg	264
STERILITY: Sterile (no growth)		
BACTERIAL SCREEN	PASSED	PASSED
FUNGAL SCREEN	PASSED	PASSED
MEA 1-cell to expanded blastocyst in 96 h	≥ 80 %	PASSED
ENDOTOXIN	< 0.25 EU/mL	PASSED

